



9079 Tyler Blvd., Mentor, OH 44060

Product name: xm18s 25

FEATURES

- Very good initial adhesion
- Gamma sterilizable
- Good for skin attachment applications
- Gentle adhesion on skin
- Good shear adhesion property
- Good for catheter fixation

COMPOSITION

Backing: 2.0 mil clear siliconized polyethylene release liner

Adhesive: 25 mil hydrocolloid adhesive

Release liner: 3.4 mil 60# SCK paper release liner

TYPICAL PROPERTIES

Test	Method	Typical Average Value
Probe Tack	Modified ASTM D2979	550-1,020 grams
Static Shear	Modified ASTM D3654A	>1,000 minutes
Saline Absorption	Xennovate Internal Method	6-20%
Liner Release	ASTM D 5357	60-120 grams

Test method details are available upon request. Values listed here are typical average values and should not be interpreted as a specification.

ROLL DESCRIPTION

Standard rolls are 9.25" wide, 450 feet long, 6" core.

RECOMMEDED STORAGE

Store in sealed, poly bags at normal room temperature (65 to 75 ° F) and dry environment. Avoid storage in temperatures above 85°F. Keep rolls in the end-cap suspended packaging supplied by Xennovate Medical[®].

Shelf Life

To obtain best performance, use this product within 12 months from date of manufacture.